

# Radius T<sup>o</sup>™

## Continuous Thermometer

en

### USER'S GUIDE

Download the Masimo Radius T<sup>o</sup> app before using this product. For help, go to [www.masimo.com/radius-t](http://www.masimo.com/radius-t) for additional tips and tutorials, a full list of supported devices, warranty, troubleshooting, and customer support.

#### INDICATIONS

Radius T<sup>o</sup> is intended for the continuous noninvasive measurement of body temperature for use on Adult and Pediatrics, 5 years of age or older, in hospitals, hospital-type facilities and home environments.

Note: Radius T<sup>o</sup> is not FDA approved or cleared.

#### CONTRAINDICATIONS

Radius T<sup>o</sup> sensors are contraindicated for people who exhibit allergic reactions to adhesive tape.

#### DESCRIPTION

Radius T<sup>o</sup> sensors are battery powered, disposable sensors that are designed to continuously provide body temperatures that are approximations of oral temperatures. The sensors are adhered to the user's skin to continuously transmit temperature measurement data via Bluetooth communication to the Masimo Radius T<sup>o</sup> application.

Note: Radius T<sup>o</sup> sensors are to only be used with compatible devices or applications. Verify compatibility before use to ensure the sensor functions properly.

#### SAFETY PRECAUTIONS

- Do not self-diagnose or self-medicate on the basis of the measurements. Always consult your doctor.
- The Masimo sensor is to only be used with a Masimo authorized app. Connection to other devices may not allow the sensor to work properly.
- Do not use the sensor if it has visible defects, discoloration, exposed wiring or appears damaged. Otherwise the sensor may not work properly.
- Avoid placing the sensor over compromised skin, excessive hair, implants, ports, skin fillers or scar tissue, as this may result in incorrect readings.
- Do not apply over or near pacemakers to avoid any potential interference from the Bluetooth communication.
- Radius T<sup>o</sup> should not be used near electrical equipment that may affect the sensor's ability from working properly.
- Check the sensor site to ensure skin integrity and to avoid damage or irritation to the skin.
- Periodically check the sensor site for proper adhesion to minimize the risk of incorrect or no readings.
- Avoid direct heating or cooling of the Radius T<sup>o</sup> sensor. Localized temperature exposure of the sensor may result in no or incorrect readings.
- Incorrect readings may be caused by sensors that are not placed on an appropriate application site.
- Radius T<sup>o</sup> may not reflect the actual body temperature when used on people undergoing treatments that may alter their normal temperature regulation (e.g. therapeutic hypothermia, antipyretics).
- Rapid or large changes in ambient temperature may cause no or incorrect readings.
- Do not modify or alter the sensor in any way. Alteration or modification may affect performance and/or accuracy.
- To prevent damage, do not soak or immerse the sensor in any liquid solution.
- Do not attempt to sterilize by irradiation, steam, autoclave, or ethylene oxide as it will damage the sensor.
- Do not attempt to reuse on multiple people. Reprocessing, reconditioning, or recycling may damage the sensor.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Radius T<sup>o</sup>, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Keep the Radius T<sup>o</sup> away from electrical equipment that emits radio frequencies to minimize radio interference. Radio interference may result in no or inaccurate readings.
- Do not use the sensor during MRI scanning or in a MRI environment as it may result in physical harm.

#### INSTRUCTIONS

- Remove sensor from pouch and pull tab to activate battery. Refer to Fig. 1.
- Pair the sensor to your phone
  - Download and open the Masimo Radius T<sup>o</sup> app.
  - Make sure your phone Bluetooth is turned on. Once battery tab is removed, your sensor is available for Bluetooth pairing. Refer to Fig. 2.
  - A solid blue light indicates the sensor is connected. (See LIGHT INDICATOR GUIDE/TROUBLESHOOTING)
  - Check your phone application to ensure the sensor is communicating correctly.

Note: Readings may take up to 15 minutes to appear on the phone application after first applying the sensor.
  - Periodically check your sensor or phone application for a solid blue light to confirm that it is connected.
- Select a site for application
  - Choose a site on the upper chest below the left collar bone. Skin should be hair-free, clean and dry before applying the sensor.

#### 4) Remove release liner from the sensor

- Remove plastic film from the back of the sensor. Refer to Fig. 3.

Note: Avoid contact with the exposed sensor adhesive.

#### 5) Apply the sensor

- Place the sensor on your selected application site. Refer to Fig. 4.

Note: Ensure that the skin is not stretched in any way and that there are no skin folds under the sensor pad.
- Press around the adhesive to ensure it is secure.

#### 6) Reapply the sensor

- Note: The sensor is designed to be removed and reapplied no more than one (1) time over the life of the product.
- Clean and dry the sensor application site.
  - Gently wipe the exposed sensor adhesive with an alcohol wipe and allow to dry to restore the adhesive properties.
  - Follow Step 5 from above to re-apply the sensor.

#### 7) Remove the sensor

- Peel gently to remove the sensor from the skin.

Note: Disposal of Product: Comply with local laws in the disposal of the sensor, battery and its accessories.

#### LIGHT INDICATOR GUIDE/TROUBLESHOOTING

Color	Sensor	Description	Next steps
No light	---	• Sensor power is off.	• Confirm battery pull tab has been removed to activate the battery. • Replace the sensor.
Green	flashing	• Sensor is on and waiting to pair with your mobile phone.	• Follow instructions to pair with your mobile phone.
Blue	flashing	• Sensor is waiting for you to confirm that desired sensor was paired to your mobile phone.	• Verify the sensor is properly attached to the wearer's body so that mobile phone can receive data.
	solid	• Successful pairing of sensor and mobile phone. • Mobile phone successfully receiving data.	
Orange	flashing	• Low sensor battery	• Consider replacing the sensor.
Red	flashing	• Depleted sensor battery • Hardware or sensor failure, sensor blinking board failure code	• Replace the sensor.

For additional help, contact Masimo Customer Support at (800) 916-1270. Local contact information can be found at [www.masimo.com/radius-t](http://www.masimo.com/radius-t).

#### PRODUCT PERFORMANCE

The Radius T<sup>o</sup> sensors have the following specifications:

Temperature measurement accuracy	±0.1°C in the range of 25°C to 43°C
Application Site	Upper Chest, below the left collarbone
Product Use/Battery Life	Minimum of 8 days, (192 hours) of continuous run time

The laboratory accuracy of Radius T<sup>o</sup> is ±0.1°C (0.18°F) for an input surface temperature range of 25°C to 43°C (77°F to 109.4°F).

Radius T<sup>o</sup> has been validated on 112 subjects, 5 years of age or older, against a reference clinical thermometer. Results have shown a clinical bias of -0.17°C (-0.30°F) with limits of agreement ≤ 0.96°C (1.73°F).

#### ENVIRONMENTAL

Storage/Transport Temperature	-20°C to 50°C @ ambient humidity
Operating Temperature	10°C to 40°C @ ambient humidity
Storage/Transport Humidity	10% RH to 95% RH (non-condensing) @ ambient temperature
Operating Humidity	10% RH to 95% RH (non-condensing) @ ambient temperature
Atmospheric Pressure	700 to 1060 hPa @ ambient temperature and humidity

# Radius T<sup>o</sup>™

## Continuous Thermometer

 Manufacturer:



Masimo Corporation  
52 Discovery  
Irvine, CA 92618  
USA  
[www.masimo.com](http://www.masimo.com)

© 2020 Masimo Corporation

301864/10890A-1020

## WIRELESS TECHNOLOGY INFORMATION

Type	Bluetooth Low Energy
Data Transmission Rate	Minimum packet rate of 0.0167 Hz (1/60 Hz)
Max. Output Power	(EIRP): 9.9 dBm
Modulation Type	GFSK
Frequency Range	2402-2480 MHz
Antenna Peak Gain	+5.67dBi

FCC ID are as follows: FCC ID: VKF-RADIUST, IC ID: 7362A- RADIUST

**CAUTION:** In order to maintain Bluetooth connectivity with the host device ensure that Radius T° is within specified distance and line of sight of the host device.

RF Radiation Exposure Statement: This equipment has been exempted from FCC RF radiation exposure testing and IC RSS 102 RF radiation exposure limits set forth for an uncontrolled environment.

**Note:** This device complies with part 15 of FCC Rules and Industry Canada's license-exempt RSSs. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Consult the dealer or an experienced radio/TV technician for help.

**Note:** When using Radius T° consideration should be taken to local government frequency allocations and technical parameters to minimize the possibility of interference to/from other wireless devices.


**Note:** Change or modifications that are not expressly approved by the manufacturer could void the user's authority to operate the equipment.

**Note:** The frequency bands of this device (2.4 GHz) are only for indoor use, in accordance with international telecommunication requirements.

## RECOMMENDED SEPARATION DISTANCES

RECOMMENDED SEPARATION DISTANCE BETWEEN PORTABLE AND MOBILE RF COMMUNICATION EQUIPMENT AND THE ME EQUIPMENT		
The ME Equipment is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ME Equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ME Equipment as recommended below, according to the maximum output power of the communication equipment.		
RATED MAXIMUM OUTPUT POWER OF TRANSMITTER (W)	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER (M)	
	80 MHz to 800 MHz $d = 1.17 \cdot \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33 \cdot \sqrt{P}$
0.01	0.12	0.23
0.1	0.37	0.74
1	1.17	2.33
10	3.7	7.37
100	11.7	23.3
For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		
<b>Note 1:</b> At 80 MHz and 800 MHz, the higher frequency range applies.		
<b>Note 2:</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS		
The ME Equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the ME Equipment should assure that it is used in such an environment.		
EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF Emissions CISPR 11	Group 1	The ME Equipment must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF Emissions CISPR 11	Class B	Suitable for use in all establishments, including domestic environments.

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY			
The ME Equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the ME Equipment should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 15 kV air	+/- 8 kV contact +/- 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50 / 60 Hz) magnetic field. IEC 61000-4-8	30 A/m	30 A/m	Guidance - Power frequency magnetic fields should be at levels characteristic of typical location in a typical hospital environment.
Portable and mobile RF communications equipment should be used no closer to any part of the ME Equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	RECOMMENDED SEPARATION DISTANCE
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	10 V/m	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the compliance level in each frequency range b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

**Note 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ME Equipment is used exceeds the applicable RF compliance level above, the ME Equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ME Equipment.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

TEST SPECIFICATIONS FOR ENCLOSURE PORT IMMUNITY TO RF WIRELESS COMMUNICATION EQUIPMENT						
TEST FREQUENCY	BAND (A) (MHZ)	SERVICE (A)	MODULATION (B)	MAXIMUM POWER (W)	DISTANCE (M)	IMMUNITY TEST LEVEL (V/M)
385	380-395	TETRA 400	Pulse modulation (b) 18 Hz	1.8	0.3	27
450	430-470	GMR5 460, FRS 460	FM (c) +/- 5 kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation (b) 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation (b) 18 Hz	2	0.3	28
870						
930						
1 720	1 700-1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 35; UMTS	Pulse modulation (b) 217 Hz	2	0.3	28
1 845						
1 970						
2 450	2 400-2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation (b) 217 Hz	2	0.3	28
5 240	5 100-5 800	WLAN 802.11 a/n	Pulse modulation (b) 217 Hz	0.2	0.3	9
5 500						
5 785						
<b>Note:</b> If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
a) For some services, only the uplink frequencies are included.						
b) The carrier shall be modulated use a 50% duty cycle square wave signal.						
c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						

The following symbols may appear on the product or product labeling:

SYMBOL	DEFINITION	SYMBOL	DEFINITION	SYMBOL	DEFINITION
	Follow instructions for use		Separate collection for electrical and electronic equipment (WEEE).		Mark of conformity to European Medical Device Directive 93/42/EEC
	Consult instructions for use		Lot code		Authorized representative in the European community
	Manufacturer		Catalogue number (model number)		Non-sterile
	Date of Manufacture YYYY-MM-DD		Masimo reference number		Not made with natural rubber latex
	Use By YYYY-MM-DD		Storage humidity Limitation		Body weight
	Do not re-use/Single patient use only		Atmospheric pressure limitation		Storage temperature range
	Federal Communications Commission (FCC) Licensing		Do not use if package is damaged and consult instructions for use		Keep dry
	Medical device		Bluetooth		Protection from ingress of particulates and water spray from any direction
	Identifies unit has been registered as a radio device		Unique device identifier		Instructions/Directions for Use/Manuals are available in electronic format @ <a href="http://www.masimo.com/radius-t">www.masimo.com/radius-t</a> <b>Note: eIFU is not available in all countries.</b>

Patents: <http://www.masimo.com/patents.htm>

Masimo, and Root are federally registered trademarks of Masimo Corporation.

Radius™ is a trademark of Masimo Corporation.

All other products, logos, or company names mentioned herein may be trademarks and/or registered trademarks of their respective companies.

